

EchoCLEAN™ RNA Kit – Protocol

for single-step removal of inhibitors and impurities from RNA solutions.

This protocol has been developed to eliminate impurities like inhibitors, salts and nucleotides from RNA solutions.

Materials and equipment needed.

- 80 to 110 µL of RNA sample for RNA samples less than 80 µL, dilute with RNase-free water or Tris buffer to a minimum of 80 µL.
- Microcentrifuge with rotor for 1.5 mL and 2 mL reaction tubes.
Important: Switch centrifuge to relative centrifugal force, rcf (x g); if this is not possible, please use formula below* to calculate the conversion of round per minute (rpm) into rcf.
- Vortexer
- One reaction tube (2 mL) per sample for column preparation
- One reaction tube (1.5 mL) per sample for elution and collection of the purified RNA.
- Pipets for 10 µL and 200 µL scale, corresponding pipet tips
- For fastest procedure (PROTOCOL 1): **Cap Puncher** (BioEcho product no. 050-001-001)

Preparation before starting.

- Set the microcentrifuge to **1,000 x g**.
- **Important:** Switch to relative centrifugal force, rcf (x g, not rpm)

Product use limitation

The EchoCLEAN RNA Kits are for research use only. They have not been registered or authorized to be used for diagnosis, prevention, or treatment of a disease.

* Most centrifuges offer the choice between rpm and g-force (rcf); if not, calculate the rpm matching the g-force using the formula: $\text{rpm} = 1,000 \times \sqrt{\frac{g}{1.12 \times r}}$, where r = radius of rotor in mm. and g the required g-force. E. g., with a radius of 150 mm, the corresponding rpm to 1,000 x g is approx. 2,400 rpm.

PROTOCOL 1: Purification using the Cap Puncher




BioEcho Cap Puncher

Column preparation

1. Vortex the **EchoCLEAN RNA Spin Column** briefly and place into a **2 mL** reaction tube. Let stand for about 5 min.
2. Use of the **Cap Puncher** (scan QR code to watch a video): Punch a hole into the column cap and lift the column together with the **Cap Puncher** out of the 2 mL collection tube. Snap off bottom closure of the column and detach the Cap Puncher by twisting while pulling out.
Place the punched spin column back into the 2 mL reaction tube.
3. Centrifuge for **1 min at 1,000 x g**. Discard the flow-through volume (void volume) collected in the 2 mL reaction tube.
4. Place the prepared spin column into a new **1.5 mL** reaction tube for elution and place back into the rack.
Continue with cleanup of RNA.

Cleanup of RNA.


5. Transfer **80 – 110 µL of the RNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below:



Insert pipet tip vertically through the hole in the column cap and pipet the sample slowly (~5 sec) into the column.

Note:

- Do not punch the resin bed with the pipette tip during loading of lysate!
6. Centrifuge **1 min at 1,000 x g**. The purified RNA elutes into the 1.5 mL elution tube and can be immediately applied in downstream applications.

The eluted RNA can be used immediately or stored at 4 °C or for long-term storage at – 80 °C. For spectrophotometric analysis, use **Low-TE Buffer**  as blank.


PROTOCOL 2: Purification **without** the use of a Cap Puncher

Column preparation

1. Vortex the **EchoCLEAN RNA Spin Column** briefly and place into a **2 mL** reaction tube. Let stand for 5 min.
2. Loosen the screw cap of the spin column **half a turn** and **snap off the bottom closure**.
Important: Do not close the screw cap of the spin column. The screw cap must stay loosened **half a turn** to avoid vacuum generation. Place the column back into the 2 mL collection tube and both into the centrifuge.
3. Centrifuge for **1 min at 1,000 x g**. Discard the flow-through volume (void volume) collected in the 2 mL reaction tube.
4. Place the prepared spin column into a new **1.5 mL** reaction tube for elution of the sample RNA and place back into the rack. Continue with cleanup of RNA.

Cleanup of RNA.

5. Transfer **80 – 110 µL of the RNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below.




Open cap and pipet the sample slowly (~5 sec) onto the center of the resin bed of the prepared spin column. Close screw cap and loosen again half a turn.

Important: Do not close the screw cap of the spin column tightly!

Note:

- Do not punch the resin bed with the pipette tip during loading of lysate!
6. Centrifuge **1 min at 1,000 x g**. The purified RNA elutes into the 1.5 mL elution tube and can be immediately applied in downstream applications.

The eluted RNA can be used immediately or stored at -80°C . For spectrophotometric analysis, use **Low-TE Buffer**  supplied with the kit as blank.

EchoCLEAN™ RNA Kit

**for single-step removal of AND
impurities from RNA solutions**

| | | |
|--------------|--|-------------------------|
| Product no. | 020-002-050-010 (10) | 020-002-050-050 (50) |
| Kit contents | EchoCLEAN RNA Spin Columns and Low-TE Buffer | |


Quick PROTOCOL (please read protocol first)

Column preparation.

- Vortex **EchoCLEAN spin column** and place in a **2 mL** tube. Let stand for 10 min.
- Punch a hole in the cap with the Cap Puncher, and break off bottom closure (scan QR code to watch a video).
- Place spin column back into 2 mL tube.
- Centrifuge **1 min at 1,000 x g** to elute column buffer.
- Place column in a **1.5 mL** tube.

Cleanup of RNA.

- Transfer **RNA sample (max. 110 µL)** by pipetting **slowly** through cap hole (scan QR code to watch a video)
- Centrifuge **1 min at 1,000 x g** to elute RNA into **elution tube**.
- Eluted RNA is ready to use.



BioEcho Life Sciences GmbH
Nattermannallee 1
50829 Köln/Cologne, Germany

Phone: +49 221 998897-0
E-Mail: contact@bioecho.de
www.bioecho.com