

EchoCLEAN Organic Solvent DNA CleanUp Kit- Protocol

for single-step removal of inhibitors and a wide range of impurities from DNA solutions.

This protocol has been developed to remove impurities like inhibitors, salts, nucleotides, TRIzol™, phenol, chloroform and other organic solvents from DNA solutions.

Materials and equipment needed

- 80 to 110 µl of DNA sample. For DNA samples less than 80 µl, dilute with Tris buffer to a minimum of 80 µl
- Microcentrifuge with rotor for 1.5 and 2 ml reaction tubes
Important: Switch centrifuge to *relative centrifugal force, rcf* (x g*); if this is not possible please use formula below* to calculate the conversion of round per minute (rpm) into rcf.
- Vortexer
- One reaction tube (2 ml) per sample for column preparation
- One reaction tube (1.5 ml) per sample for elution and collection of the purified DNA
- Pipets for 10 µl and 200 µl scale, corresponding pipet tips
- For fastest procedure (PROTOCOL 1): Cap Puncher (BioEcho product no. 050-001-001)

Preparation before starting

- Set the microcentrifuge to **1000 x g***
- **Important: Switch to relative centrifugal force, rcf (x g*, not rpm)**

Product use limitation

The EchoCLEAN Organic Solvent CleanUp Kits are for research use only. They have not been registered or authorized to be used for diagnosis, prevention or treatment of a disease.

* Most centrifuges offer the choice between rpm and g-force (rcf); if not, calculate the rpm matching the g-force using the formula: $\text{rpm} = 1,000 \times \sqrt{\left(\frac{g}{1.12 \times r}\right)}$, where r = radius of rotor in mm. and g the required g-force.

E. g., with a radius of 150 mm, the corresponding rpm to 1,000 x g is approx. 2,400 rpm.

PROTOCOL 1: Purification using the Cap Puncher



BioEcho Cap Puncher

Column preparation

1. Vortex the **EchoCLEAN Organic Solvent Spin Column** briefly and place into a **2 ml** reaction tube.
Let stand for about 5 min.
2. Use of the Cap Puncher (scan QR code to watch a video): Punch a hole into the column cap and lift the column together with the Cap Puncher out of the 2 ml collection tube. Snap off bottom closure of the column and detach the Cap Puncher by twisting while pulling out.
Place the punched spin column back into the 2 ml reaction tube
3. Centrifuge for **1 min at 1000 x g***. Discard the flow-through volume (“void volume”) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new **1.5 ml** reaction tube for elution of the purified DNA and place back into the rack.
Continue with “Clean up of DNA”.

Clean up of DNA


5. Transfer **80 – 110 µl of the DNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below:



Insert pipet tip vertically through the hole in the column cap and pipet the sample slowly (~5 sec) into the column.

Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge **1 min at 1000 x g***. The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications.

The eluted DNA can be used immediately or stored at 4°C or for long-term storage at -20°C. For spectrophotometric analysis, use 1x Tris Buffer  supplied with the kit as blank.

EchoCLEAN Organic Solvent DNA CleanUp Kit

for 1-step removal of inhibitors, TRIzol™, organic solvents and a wide range of impurities from DNA solutions


PROTOCOL 2: Purification **without** a Cap Puncher

Column preparation

1. Vortex the **EchoCLEAN Organic Solvent Spin Spin Column** briefly and place into a **2 ml** reaction tube. Let stand for 5 min.
2. **Loosen** the screw cap of the spin column **half a turn** and **snap off the bottom closure**.
Important: Do not close the screw cap of the spin column. The screw cap must stay loosened **half a turn** to avoid generation of a vacuum. Place the column back into the 2 ml collection tube and both into the centrifuge.
3. Centrifuge for **1 min at 1000 x g***. Discard the flow-through volume (“void volume”) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new **1.5 ml** reaction tube for elution of the sample DNA and place back into the rack. Continue with “Clean up of DNA.”

Clean up of DNA

5. Transfer **80 – 110 µl of the DNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below



Open cap and pipet the sample slowly (~5 sec) onto the center of the resin bed of the prepared spin column. Close screw cap and loosen again half a turn.
Important: Do not close the screw cap of the spin column tightly!

Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge **1 min at 1000 x g***. The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications

The eluted DNA can be used immediately or stored at 4°C or for long-term storage at -20°C. For spectrophotometric analysis, use 1x Tris Buffer **T** supplied with the kit as blank.

Product no. (rxn's)	020-002-040-010 (10)	020-002-040-010 (50)	020-002-040-250 (250)
Kit contents	EchoCLEAN Organic Solvent DNA CleanUp Spin Columns, 1x Tris Buffer		


Quick PROTOCOL (please read protocol first)

Column preparation

- Vortex **EchoCLEAN spin column** and place in a **2 ml** tube
Let stand for 10 min
- Punch a hole in the cap with the Cap Puncher, and break off bottom closure (scan QR code to watch a video)
- Place spin column back into 2 ml tube
- Centrifuge **1 min at 1000 x g*** to elute column buffer
- Place column in a **1.5 ml** tube

Clean up of DNA

- Transfer **DNA sample (max. 110 µl)** by pipetting **slowly** through cap hole (scan QR code to watch a video)
- Centrifuge **1 min at 1,000 x g*** to elute DNA into **Elution tube**
- Eluted DNA is ready to use



BioEcho Life Sciences GmbH
Nattermannallee 1
50829 Köln (Cologne)/Germany

Phone: +49 (0) 221-99 88 97-0
E-Mail: contact@bioecho.de
www.bioecho.de