

# EchoCLEAN Organic Solvent DNA CleanUp Kit- Protocol

for single-step removal of inhibitors and a wide range of impurities from DNA solutions.

This protocol has been developed to remove impurities like inhibitors, salts, nucleotides, TRIzol™, phenol, chloroform, and other organic solvents from DNA solutions.

## Materials and equipment needed

- 80 to 110 µl of DNA sample. For DNA samples less than 80 µl, dilute with Low-TE Buffer to a minimum of 80 µl.
- Microcentrifuge with rotor for 1.5 ml and 2 ml reaction tubes.  
**Important:** Switch centrifuge to *relative centrifugal force, rcf* (x g\*); if this is not possible, please use formula below\* to calculate the conversion of round per minute (rpm) into *rcf*.
- Vortex.
- One reaction tube (2 ml) per sample for column preparation.
- One reaction tube (1.5 ml) per sample for elution and collection of the purified DNA
- Pipets for 10 µl and 200 µl scale, corresponding pipet tips
- For fastest procedure (PROTOCOL 1): Cap Puncher (BioEcho product no. 050-001-001)

## Preparation before starting

- Set the microcentrifuge to **1000 x g\***
- **Important: Switch to relative centrifugal force, rcf (x g\*, not rpm)**

## Product use limitation

The EchoCLEAN Organic Solvent CleanUp Kits are for research use only. They have not been registered or authorized to be used for diagnosis, prevention, or treatment of a disease.

\* Most centrifuges offer the choice between rpm and g-force (rcf); if not, calculate the rpm matching the g-force using the formula:  $\text{rpm} = 1,000 \times \sqrt{\left(\frac{g}{1.12 \times r}\right)}$ , where r = radius of rotor in mm. and g the required g-force.

E. g., with a radius of 150 mm, the corresponding rpm to 1,000 x g is approx. 2,400 rpm.

## PROTOCOL 1: Purification using the Cap Puncher



BioEcho Cap Puncher


### Column preparation

1. Vortex the **EchoCLEAN Organic Solvent Spin Column** briefly and place into a **2 ml** reaction tube.  
Let stand for about 5 min.
2. Use of the Cap Puncher (scan QR code to watch a video): Punch a hole into the column cap and lift the column together with the Cap Puncher out of the 2 ml collection tube. Snap off bottom closure of the column and detach the Cap Puncher by twisting while pulling out.  
Place the punched spin column back into the 2 ml reaction tube
3. Centrifuge for **1 min at 1000 x g\***. Discard the flow-through volume (void volume) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new **1.5 ml** reaction tube for elution of the purified DNA and place back into the rack.

**Continue with DNA Cleanup.**

### DNA Cleanup


5. Transfer **80 – 110 µl of the DNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below:



Insert pipet tip vertically through the hole in the column cap and pipet the sample slowly (~5 sec) into the column.

### Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge **1 min at 1000 x g\***. The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications.

The eluted DNA can be used immediately or stored at 4°C or for long-term storage at -20°C. For spectrophotometric analysis, use Low-TE Buffer  supplied with the kit as blank.

# EchoCLEAN Organic Solvent DNA CleanUp Kit

for 1-step removal of inhibitors, TRIzol™, organic solvents and a wide range of impurities from DNA solutions


## PROTOCOL 2: Purification **without** a Cap Puncher

### Column preparation

1. Vortex the **EchoCLEAN Organic Solvent Spin Column** briefly and place into a **2 ml** reaction tube. Let stand for 5 min.
2. **Loosen** the screw cap of the spin column **half a turn** and **snap off the bottom closure**.  
**Important: Do not close the screw cap of the spin column.** The screw cap must stay loosened **half a turn** to avoid generation of a vacuum. Place the column back into the 2 ml collection tube and both into the centrifuge.
3. Centrifuge for **1 min at 1000 x g\***. Discard the flow-through volume (void volume) collected in the 2 ml reaction tube.
4. Place the prepared spin column into a new **1.5 ml** reaction tube for elution of the sample DNA and place back into the rack.  
**Continue with DNA Cleanup.**

### DNA Cleanup


5. Transfer **80 – 110 µl of the DNA sample** to the prepared **EchoCLEAN Spin Column** from step 4 as illustrated below



Open cap and pipet the sample slowly (~5 sec) onto the center of the resin bed of the prepared spin column. Close screw cap and loosen again half a turn.  
**Important: Do not close the screw cap of the spin column tightly!**

#### Note:

- Do not punch pipette tip into the resin bed during loading of lysate!
6. Centrifuge **1 min at 1000 x g\***. The purified genomic DNA elutes into the 1.5 ml elution tube and can be immediately applied in downstream applications

The eluted DNA can be used immediately or stored at 4°C or for long-term storage at -20°C. For spectrophotometric analysis, use 1x Tris Buffer  supplied with the kit as blank.

|                     |   |                         |                          |
|---------------------|---|-------------------------|--------------------------|
| Product no. (rxn's) | 020-002-040-010<br>(10)   | 020-002-040-050<br>(50) | 020-002-040-250<br>(250) |
| Kit contents        | EchoCLEAN Organic Solvent DNA CleanUp Spin Columns, and Low-TE Buffer |                         |                          |


### Quick PROTOCOL (please read protocol first)

#### Column preparation

- Vortex **EchoCLEAN spin column** and place in a **2 ml** tube  
Let stand for 10 min.
- Punch a hole in the cap with the Cap Puncher, and break off bottom closure (scan QR code to watch a video).
- Place spin column back into 2 ml tube.
- Centrifuge **1 min at 1000 x g\*** to elute column buffer.
- Place column in a **1.5 ml** tube.

#### DNA Cleanup

- Transfer **DNA sample (max. 110 µl)** by pipetting **slowly** through cap hole (scan QR code to watch a video).
- Centrifuge **1 min at 1,000 x g\*** to elute DNA into **Elution tube**.
- Eluted DNA is ready to use.



BioECHO Life Sciences GmbH  
Nattermannallee 1  
50829 Köln (Cologne)/Germany

Phone: +49 (0) 221-99 88 97-0  
E-Mail: [contact@bioecho.de](mailto:contact@bioecho.de)  
[www.bioecho.de](http://www.bioecho.de)